

FECAL BACTERIAL CULTURE 1 (F -BaktVi1)

F -BaktVi1 includes F -SalmVi, F -ShigVi, F -YersVi and F -CampVi

- Collect some stool into a clean bowl.
- Use the spoon on the lid of the sample container to collect a small amount (one teaspoon) of sample. Close the container carefully.
- Open the transport tube pack. Remove the white cap from the sample tube. Use the sampling swab and collect stool sample especially from an area of blood and mucus. Put the sampling swab to bottom of the transportation tube under the gel. Do not break or cut the swab. The blue cap at the end of the swab will seal the tube.
- Write down your name, social security number and the date of collection on stickers. Attach the stickers onto the side of the sample tube.
- Put the sample container and tube in the plastic bag and seal the bag carefully.
- Carefully fill out all the fields in this form. Deliver the samples and this form back to the laboratory.

Please fill in carefully following fields:

Name _____

Social security number / date of birth: _____

Date of collection: ____/____/____

Have you travelled abroad within the past 6 months? Yes No

Country and time of travel: _____

Have you experienced abdominal symptoms during or after the travel? Yes No

Inquiries: