



Fecal blood test (F -hHb)




1. Collect the stool on clean toilet paper or in a clean container (e.g., a disposable plate or cup). Do not take the sample from the toilet bowl.



2. Open the green cap of the tube. Do not touch the clear cap. The sampling stick is attached to the cap.



3. . Collect the stool sample by inserting the stick into four different spots. Finally, run the stick through the stool from end to end in a crosswise motion.



4. The sample is sufficient when the grooves on the stick are covered with stool.



5. Place the green cap back on the tube and close it securely. Mix the contents by turning the tube several times. Do not open the tube after this.



6. Place the sample tube in the sealable bag and close the bag securely. Deliver the sample and this filled form back to the laboratory.

Please fill in the following fields carefully:

Name: _____

Personal ID: _____

Date of sample collection: ___/___/___ Time:__:__

Note! Do not touch the liquid inside the tube, and do not open the clear cap of the tube. Avoid spilling the contents of the tube on skin, mouth or eyes. However, if this happens, immediately rinse the area thoroughly. Do not drink the liquid. Keep out of reach of children..