

Calprotectin (F -Calpro)

- Collect some stool into a clean bowl.
- Use the spoon on the lid of the sample container to collect the sample (teaspoonful, no more than half the volume of the sample container). Close the container carefully.
- Write down your name, personal identity number and the date of collection on a label sticker. Attach the sticker onto the side of the container. Place the container into a sample bag and close the bag. Carefully fill out all the fields in this form.
- Deliver the sample and this completed form to the laboratory where you received the sampling supplies.

Please fill in carefully the following fields:

Name: _____

Social security number/date of birth: _____

Date of collection: ____/____/____

Inquiries: